




Low Cross Talk 384 Well Polypropylene Plate



- Low cross talk and minimal binding
- Ideal for luminescence applications
- Highly reflective surface
- Low autophosphorescence
- Minimum capillary action
- Standard height and foot print

Nunc Low Cross Talk 384 Well Polypropylene Plate

This plate has some unique features, which make it an obvious plate to use in luminescence assays including SPA (Scintillation Proximity Assays)*.

- Developed to fill a requirement on the market for a light dense, low cross talk, highly reflective and low binding plate.
- Autophosphorescence from daylight exposure is particularly low and fast-fading due to a selected resin composition.
- Polypropylene plates are preferred to polystyrene plates for homogeneous assays (liquid phase assays) due to their lower adsorption and higher chemical resistance.

* SPA technology is covered by patents held by Nycomed Amersham plc.

- The Nunc plate exhibits low autofluorescence, a clear advantage in a variety of applications including fluorescence assays.
- With the same dimensions as the other Nunc 384 Well polypropylene plates and a total volume of 120 μ l this plate conforms to the recommended standard.
- In addition the Nunc Low Cross Talk Plate has improved stacking alignment. It is compatible with automated equipment e.g. stackers, workstations and liquid dispensers.
- The inherent hydrophobic properties and rounded square wells of the Nunc polypropylene plates minimise capillary action (wicking), which is particularly a problem in plates with a very hydrophilic surface (fig. 2).

Cross-talk comparison

The Nunc Low Cross Talk 384 Well Polypropylene Plate is compared with a plate of a competitive brand (fig. 1). As can be seen in the graph below, the Nunc Low Cross Talk Plate shows lower cross talk and a higher signal than the competitor plate.




Fig. 1
On the left axis the percentage cross talk from adjacent wells is shown as a proportion of the total signal on the right axis in IOD (Integrated Optical Density) units measured by SPA in a LEADseeker™ instrument.

Capillary action comparison




Fig. 2.
Wells filled with 100 μ l aqueous colour solution show minimal capillary action (wicking effect) in a Nunc Low Cross Talk 384 Well Polypropylene Plate (to the left) in contrast to the effect in a polystyrene plate (to the right). Capillary action increases the risk of overflow, particularly in conjunction with sealing.

Ordering information:

Nunc Low Cross Talk 384 Well Polypropylene Plate

Cat. No.	264675
Colour	White
Sterilized	–
Units per pack/case	20/120