


## MVE Lab Series


The Lab Series cryogenic liquid dewars are named for their acceptance in laboratories and medical facilities worldwide. These high-efficiency, super insulated dewars are the most convenient, economical way to store and dispense liquid nitrogen. Many lab units can be fitted with pouring spouts, pressurized dispensing devices or dippers to aid in the transfer of liquid nitrogen.

### Features Include:

- Designed for efficient storage of liquid nitrogen
- Low liquid nitrogen consumption
- Convenient lightweight package



## Tank Features



	LAB 4	LAB 5	LAB 10	LAB 20	LAB 30	LAB 47	LAB 50
<b>Performance</b>							
<b>LN2 Capacity w/o Inventory</b> L est.	4	5	10	20	32	47	50
<b>Static Evaporation Rate*</b> L/day	0.19	0.15	0.18	0.18	0.22	0.39	0.49
<b>Unit Dimensions</b>							
<b>Neck Opening</b> in. (mm)	1.40 (35)	2.18 (56)	2.18 (56)	2.18 (56)	2.50 (64)	5.00 (127)	2.50 (64)
<b>Overall Height</b> in. (mm)	16.8 (426)	18.2 (462)	21.5 (546)	24.5 (622)	24.0 (610)	26.5 (673)	30.5 (775)
<b>Usable Height</b> in. (mm)	7.8 (198)	10.5 (266)	13.5 (343)	13.7 (348)	14.9 (378)	16.8 ( 427)	22.0 (559)
<b>Outer Diameter</b> in. (mm)	7.3 (185)	8.8 (222)	10.3 (260)	14.5 (368)	17.0 (432)	20.0 (508)	17.0 (432)
<b>Internal Diameter</b> in. (mm)	5.5 (139)	6.5 (165)	8.3 (210)	11.4 (289)	14.0 (356)	16.5 (419)	14.0 (356)
<b>Weight Empty</b> lb. (kg)	6 (2.7)	8 (4)	12 (5.4)	19 (9)	25 (11.4)	42 (19)	31 (14)
<b>Weight Full</b> lb. (kg)	13 (6)	17 (8)	31 (14)	55 (25)	82 (37.2)	120.4 (54.6)	120 (54.4)

\* Static evaporation rate and static holding time are nominal. Actual rate and holding time will be affected by the nature of container use, atmospheric conditions, and manufacturing tolerances.

Note: The Lab Series is not intended for sample storage.

**TWO Year Parts Warranty • FIVE Year Vacuum Warranty**

Conforms to MDD 93/42/EEC, the Medical Device Directive for the EU.



**Labo Baza**

**nowoczesne wyposażenie laboratorium**

ul. Topolowa 5  
62-002 Jelonek k/Poznania  
tel./fax: 61 812 57 45  
e-mail: biuro@labobaza.pl  
**www.labobaza.pl**