

Ultrapure water made easy!

The Thermo Scientific™
Barnstead™ GenPure™ Pro
water purification system
is ideal for the most
demanding laboratory
applications. This system
delivers ultrapure water from
a small footprint and can be
used on the bench or wall.



innovation.

▶ UV intensity and feed water monitoring

Continuous monitoring of UV intensity reduces the risk of false TOC values, while feed water monitoring alerts you to fluctuations in feed water quality.

flexibility.

▶ Flexible dispensing

New flexible dispenser offers a radius of 24 in (60 cm) from the system for filling larger vessels or glassware washing.

convenience.

▶ Simple self-maintenance

Aquastop quick-connect capability enables cartridge replacement within seconds—even mid-operation.

**Special savings
available!**

**Contact your sales
representative for details.**

Cost-effective ordering— start-up consumables included!

All consumables needed for start-up come standard with system. Includes feed water pressure regulator, ultrapure polishing cartridge, wall mounting bracket, sterile filter, UV lamp, ultrafilter (on applicable units) and TOC monitor (on applicable units).

Thermo
SCIENTIFIC

Special savings available!


Contact your sales representative for details.



Quick connect cartridge replacement



Exact dosing with easy push-button dispensing



Intuitive display shows important system parameters

Thermo Scientific Barnstead GenPure Pro water purification system

Description	Application	Promotional Part Number
GenPure Pro UV	Analytical	50131952PR
GenPure Pro UV TOC	Analytical	50131948PR
GenPure Pro UV/UF	Life Science	50131950PR
GenPure Pro UV/UF TOC	Life Science	50131922PR

Analytical models are ideal for applications such as inorganic and organic trace analysis, HPLC, ICP-MS, IC, and TOC Analysis.

Life science models are ideal for applications such as PCR, DNA, RNA, monoclonal antibodies, or cell culture media prep requiring nuclease removal.

TOC models provide visual read-out of TOC levels.

All consumables needed for start-up come standard with system. Includes feed water pressure regulator, ultrapure polishing cartridge, wall mounting bracket, sterile filter, UV lamp, ultrafilter (on applicable units) and TOC monitor (on applicable units).

© 2013 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Valid with participating distributors in Czech Republic, Greece, Hungary, Israel, Poland, Saudi Arabia, South Africa, Turkey and United Arab Emirates while supplies last. All promotions valid for orders placed between May 15, 2013 and December 31, 2013. "Savings" or "Discount" indicated in any promotion herein represent the maximum percentage difference between the aggregate price for individual purchase of each item contained in the promotional package at the current standard local MSRP and the local MSRP for the promotional price for each promotional package. Valid only in connection with purchase of promotional packages identified herein using the part numbers listed. All pricing information and values based on current list prices. All offers subject to Thermo Fisher standard terms and conditions of sale, and Thermo Fisher makes no additional warranties, express or implied, in connection with promotional products. Promotional products will be shipped with purchased products. Not valid in connection with any other offer or discount. Void where prohibited by law or company policy. See representative for additional information.

act now!

Contact your sales representative to order your Barnstead GenPure Pro water purification system package today!

Thermo
SCIENTIFIC
Part of Thermo Fisher Scientific