



## Homogenizer Safe-Seal

**Fisherbrand®**

Complete homogeniser with PTFE pestle, stainless steel handle and borosilicate glass vessel with PTFE stopper as splash guard.

Volume ml	Ø pestell mm	Vessel length mm	PU pcs.	Cat. No.
5	10	65	1	<b>FB 56763</b>
10	13	75	1	<b>FB 56765</b>
15	15	80	1	<b>FB 56767</b>
30	19	105	1	<b>FB 56769</b>




## Mini Homogenizer, Borosilicate Glass

**Fisherbrand®**

High quality borosilicate glass. Mini homogeniser is ideal for working with micro tissue. Delivered as complete apparatus suitable for autoclaves. Vessel made of transparent precision glass, pestle with ground head.

Volume ml	Ø pestell mm	Vessel length mm	PU pcs.	Cat. No.
0.1	3.65	31	1	<b>FB 56673</b>
1.0	5.0	55	1	<b>FB 56675</b>
3.0	9.0	60	1	<b>FB 56679</b>



## Homogenizers Potter-Elvehjem Type, Glass

**KIMBLE CHASE**  
Life Science and Research Products

Suitable for use with motor.

Type	Type	Volume ml	PU pcs.	Cat. No.
Vessel	19	1	1	<b>312 20 01</b>
Vessel	21	5	1	<b>312 21 30</b>
Vessel	22	8	1	<b>312 21 31</b>
Vessel	23	17	1	<b>312 21 32</b>
Vessel	24	45	1	<b>312 21 33</b>
Glass pestle	19	1	1	<b>312 21 34</b>
Glass pestle	21	5	1	<b>312 21 35</b>
Glass pestle	22	8	1	<b>312 21 36</b>
Glass pestle	23	17	1	<b>312 21 37</b>
Glass pestle	24	45	1	<b>312 21 38</b>



## Homogenizers Potter-Elvehjem Type, Glass Vessel and PTFE Pestle

**KIMBLE CHASE**  
Life Science and Research Products

Suitable for use with motor.

Type	Type	Volume ml	PU pcs.	Cat. No.
Vessel	18	0.5	1	<b>312 20 02</b>
Vessel	19	1	1	<b>312 21 39</b>
Vessel	20	3	1	<b>312 21 46</b>
Vessel	21	5	1	<b>312 21 40</b>
Vessel	22	8	1	<b>312 21 41</b>
Vessel	23	17	1	<b>312 21 42</b>
Vessel	24	45	1	<b>312 21 43</b>
PTFE pestle	18	0.5	1	<b>312 21 44</b>
PTFE pestle	19	1	1	<b>312 21 45</b>
PTFE pestle	20	3	1	<b>312 21 47</b>
PTFE pestle	21	5	1	<b>312 21 48</b>
PTFE pestle	22	8	1	<b>312 21 49</b>
PTFE pestle	23	17	1	<b>312 21 50</b>
PTFE pestle	24	45	1	<b>312 21 51</b>